




Automation

SeqCap EZ HyperCap Workflow enables easy automation

Vendor	Instrument	KAPA Library Prep Kit	Hyper Prep Kit	Hyper Plus Kit	SeqCap EZ v5.1	HyperCap v1.0
Perkin Elmer	Sciclone	✓	✓	✓	✓	✓
	Zephyr	✓				
	NGS Express	✓	✓		✓	
Beckmann	FX	✓	✓	✓	✓	✗
Tecan	EVO NGS	✓				
Hamilton	STAR NGS		✓	✓	✓	✗
Agilent	Bravo B	✓	✓	✓	✓	✗
Eppendorf	epMoption TMX 5075	✓	✓	✓		
 Refers to methods in development						

Automation platforms available for use with SeqCap EZ enrichment products.
Note: The availability of "Standard Solutions" is constantly changing. Contact your instrument providers or see KAPA's Automation webpage to obtain the most up-to-date list of available applications.