

Mastermix Selection Guide

Welcher Mastermix ist die richtige Wahl für welches LightCycler® System?

		LightCycler® 96 System	LightCycler® 480 System	LightCycler® PRO System	
Absolute Quantification/ Gene Expression (Relative Quantification, Two-Step RT-PCR)	TaqMan® Probes	●		●	FastStart Essential DNA Probes Master
			●		LightCycler® 480 Probes Master
		●	●	●	LightCycler® Multiplex DNA Master
	SYBR Green	●		●	FastStart Essential DNA Green Master
			●		LightCycler® 480 SYBR Green I Master
	HybProbe Probes			●	LightCycler® FastStart DNA Master HybProbe
Gene Expression (Rel. Quantification, One-Step RT-PCR)	TaqMan® Probes	●	●	●	LightCycler® Multiplex RNA Virus Master
	TaqMan® Probes	●		●	FastStart Essential DNA Probes Master
			●		LightCycler® 480 Probes Master
	HybProbe Probes			●	LightCycler® FastStart DNA Master HybProbe
			●	●	LightCycler® 480 Genotyping Master
	High Resolution Melting (HRM) Dye	●	●		LightCycler® 480 High Resolution Melting Master

Bestellinformation

Produkt	Inhalt	Bestellnummer
LightCycler® 480 SYBR Green I Master	500 Reaktionen je 20 µl	04 707 516 001
	5.000 Reaktionen je 20 µl	04 887 352 001
LightCycler® 480 Probes Master	500 Reaktionen je 20 µl	04 707 494 001
	5.000 Reaktionen je 20 µl	04 887 301 001
LightCycler® Multiplex DNA Master	200 Reaktionen je 20 µl	07 339 585 001
	1.000 Reaktionen je 20 µl	07 339 577 001
DNA Process Control Kit	600 Reaktionen je 20 µl	07 339 542 001
LightCycler® 480 Genotyping Master	384 Reaktionen je 20 µl	04 707 524 001
LightCycler® 480 High Resolution Melting Master	500 Reaktionen je 20 µl	04 909 631 001
RNA Process Control Kit	200 Reaktionen je 20 µl	07 099 622 001
	600 Reaktionen je 20 µl	07 099 592 001
LightCycler® Multiplex RNA Virus Master	200 Reaktionen je 20 µl	06 754 155 001
	1.000 Reaktionen je 20 µl	07 339 577 001
Fast Start Essential DNA Green Master	500 Reaktionen je 20 µl	06 402 712 001
	5.000 Reaktionen je 20 µl	06 924 204 001
Fast Start Essential DNA Probes Master	500 Reaktionen je 20 µl	06 402 682 001
	5.000 Reaktionen je 20 µl	06 924 492 001

Unsere Systeme



LightCycler® PRO System



LightCycler® 480 System



LightCycler® 96 System

Roche Diagnostics Deutschland GmbH
Sandhofer Straße 116
68305 Mannheim

LIGHTCYCLER und TAQMAN sind Marken von Roche.
© 2025 Roche Diagnostics. Alle Rechte vorbehalten.

www.roche.de

Sie wollen mehr erfahren?

Kontaktieren Sie uns unter:
mannheim.vertrieb-molecular@roche.com

roche.de/lightcycler

