

cobas® p 612 pre-analytical system

Technical specifications



System description

Basic System

cobas® p 612 pre-analytical system

08166145001

Basic Specifications

- Continuous loading of tubes
- Simultaneous loading of different tube types
- Prioritization of urgent samples
- Registration of samples
- Comprehensive inspection of sample quality
 - Tube type identification via tube dimensions & cap color
(used to cross-check the tube type with the test request and/or specimen)
 - Liquid level check and sample volume calculation through 3 barcode label
 - Spin status detection
 - Sample clot detection (during aliquoting)
- Decapping of tubes
- Simultaneous decapping of hemogard, screw and rubber stopper caps
- Removal of caps with a twist and pull motion to prevent build-up of aerosols
- Aliquoting according configurable criteria, including utilizing sample quality check information (up to 28 aliquots from one primary tube)
- Aliquot labelling
- Safety cross check of the aliquot barcode with the primary barcode
- Definable pipetting volume (disposable tip is not exchanged in case of multiple sample uptake)
- Sorting according configurable sorting criteria, including utilizing sample quality check information for error handling
- Flexible and freely definable input and output areas:
 - Into and out of racks
 - Out of centrifuge buckets
- Definable output areas for tubes with specific errors (e.g. BC not readable, no test requests)
- First pass and recursive workflow
- Automatic barcode alignment of tubes for reading and sorting
- Archiving of samples

Options	Possibility to connect to cobas® connection modules (CCM) with different rack configurations.	Multiple material numbers
	Bulk Loader Input Module for cobas® p 512 and cobas p 612 Bulk Loader Input Module for cobas® p 471 and cobas® p 671 cobas p 471 centrifuge unit cobas p 671 centrifuge unit	09469761001 09469460001 07962606001 07962622001
Optional Specifications	<ul style="list-style-type: none"> • Bulk loading of samples • Centrifugation of samples (connection to single or double centrifuge) 	<ul style="list-style-type: none"> • Capture of sample image • Recognition of sample quality into 5 different user-defined categories (serum ok, light hemolytic, hemolytic, lipemic, icteric)
Specifications		
Throughput	Up to 1400 tubes per hour (basic system without aliquoting)*	
Sample identification	<ul style="list-style-type: none"> • Positive sample identification via barcode • Barcode types supported: Codabar, 2 of 5 interleaved (ITF), Code 39, Code 128 (incl. ISBT 128), others** 	
Input / Output Capacity	<ul style="list-style-type: none"> • Input: 600 samples, 4 drawers, up to 30 targets • Output: 1200 samples; 8 drawers, up to 41 targets 	
Tube types	<p>Supports established tubes for hematology, coagulation, serum, plasma and urine manufactured by BD, Sarstedt, Kabe, Greiner, Terumo and others tube vendors</p> <ul style="list-style-type: none"> • Accepts most 3, 5, 7 and 10 mL plastic tubes • Uses qualified tubes (list with qualified tubes is available from your distributor) • Tubes can be qualified on request 	
Sample carriers	<p>Supports racks supplied from Roche and other vendors</p> <ul style="list-style-type: none"> • Standard racks 	<ul style="list-style-type: none"> • Analyzer racks • Archive racks • Centrifuge buckets
Aliquot tips	<ul style="list-style-type: none"> • Disposable tips for carry-over free • Standard pipetting tips, non-conductible 1000 µL (10 x 96) 04976274001 	<ul style="list-style-type: none"> • 384 tips on board in 4 tip racks • Dualfilter pipetting tips, non-conductible 1000 µL
Aliquot tubes	<ul style="list-style-type: none"> • Round bottom tube 13 x 75mm 	<ul style="list-style-type: none"> • False bottom tube 13 x 75 mm
Dimensions	<ul style="list-style-type: none"> • Width: approx. 3.14 m (123.6") • Height: approx. 1.84 m (72.4") <p>Control unit is incorporated on the system</p>	<ul style="list-style-type: none"> • Depth: approx. 1.41 m (55.5") • Weight: approx. 878 kg (1936 lbs)
Operating conditions	<ul style="list-style-type: none"> • Designed for indoor operation up to 2000 meters above sea level • Ambient temperature +15 °C to +30 °C (+59 °F to +86 °F) 	<ul style="list-style-type: none"> • Relative humidity 20 % to 80 % • Power consumption: approx. 800 W • Voltage: 230 V, 50/60 Hz (fuse 6.3 AH) and 115 V, 50/60 Hz (fuse 10 AH)
Compressed air	<ul style="list-style-type: none"> • Dry and oil free • Air pressure: min. 6.5 bar (94 psi) to max. 8.0 bar (116 psi) • Air consumption: approx. 94 l/min (24.8 gals/min), with recapper 	
Interface	<ul style="list-style-type: none"> • ASTM protocol • Network connection via RS 232 and Ethernet • Remote access for diagnosis and service 	
Safety marks	CE, ETL	

COBAS is a trademark of Roche.
© 2024 Roche

Published by:
Roche Diagnostics International Ltd
CH-6343 Rotkreuz
Switzerland

[diagnostics.roche.com](https://www.diagnostics.roche.com)

* The final throughput depends on different parameters as for example the aliquot volume or the number of aliquots per primary tube

** Other barcode types can be checked on request